

ottobock.

Agilium[®] Forte

Less Pain. More Life.


Quality for life



Pain Relief from Knee Osteoarthritis

Meet the Agilium Forte, a lightweight, wraparound knee brace designed to provide lasting pain relief for people with knee osteoarthritis.

With one pull, this brace uses the dynamic Y-shaped strap system to reduce pressure and relieve pain on the side of the knee that is most affected.




Break the cycle of knee OA

As knee osteoarthritis causes pain, individuals compensate their posture to avoid pressure on their knee for pain relief. This can cause issues in other parts of the body, such as the hips. Reduced movement in the knee decreases the supply of nutrients to the cartilage, which breaks down even more, increasing pain.

Indications

Indications for use to be determined by your physician.

- Moderate to severe knee osteoarthritis
- Inner or outer knee pain
- Rheumatoid knee arthritis
- Post-operative knee support



Treatment with the Agilium Forte

Regular use of a knee brace can help break the cycle of knee OA by stabilizing the knee joint and providing pain relief for one side (inner or outer) of the knee. The dynamic Y-shaped strap opens and unloads the affected side of the knee joint.

Overview

The Agilium Forte is ideal for individuals with moderate to severe knee osteoarthritis who prefer a rigid, lightweight brace that is comfortable and easy-to-use.



Lasting pain relief

A dynamic strap provides pain relief with just one pull



Lightweight design

Rigid shells for stability and open sleeve design for comfort



Easy to apply and remove

Three adjustable, numbered closure straps on the calf and thigh for quick application and excellent compression




Breathable material

Soft material with no back-of-knee bunching for comfort

Sizing

The Agilium Forte comes in three sizes. Simply measure the circumference of the thigh 6 inches above the center of the patella, and the calf 6 inches below.



Size	Thigh Circumference	Calf Circumference
XS-S	14"-18"	11"-14"
M-L	18"-23.5"	14"-18"
XL-3XL	23.5"-32"	18"-24"

Less Pain. More Life.

Osteoarthritis (OA) is one of the most common joint disorders worldwide.¹ As the aging population dramatically increases, more people are facing mobility challenges. When it affects the knee, osteoarthritis is one of the most frequent causes of difficulty with walking and climbing stairs.

A knee brace may help in the treatment of osteoarthritis to reduce pain, improve walking, and increase stability in patients with knee OA. Drawing from 100 years of experience in quality medical devices, Ottobock developed the Agilium Forte, a rigid knee osteoarthritis brace for everyday use, recreational activities, and sports. This addition to the Agilium portfolio of knee braces dynamically relieves knee pain with a single strap so you can feel less pain, and experience more life.

One Pull, and Go!

Get back to the activities you love with the pain alleviation you need to lead a healthy, active life. Contact your physician to determine if the Agilium Forte is the right fit for your needs. Upon initial fitting, your physician can adjust the sizing and customize the fit of the brace for maximum pain relief.



¹Grazio 2009

